



## TEKNOMAR BRAND

## PASSBOX DATASHEET

Teknomar brand passbox is manufactured to prevent unnecessary human traffic and provide fast transportation of products.

GENERAL SPECIFICATIONS		
	PASSIVE (STATIC)	ACTIVE (DYNAMIC)
HEPA FILTER	NO	YES
UV	OPTIONAL	YES
AIR CYCLING	NO	YES
VENTILATOR	NO	YES
INTERLOCK	YES	YES
GLASS TYPE	TEMPERED	TEMPERED
DOOR TYPE	DOUBLE	DOUBLE
INSIDE EQUIPM.	304 L STAINLESS STEEL	304 L STAINLESS STEEL
BODY EQUIPM.	304 L / ALUMINYUM (OPTIONAL)	304 L STAINLESS STEEL
INSIDE DIMENSIONS	60x60x60cm	60x60x60cm
CONTROL	ELECTRONIC	ELECTRONIC
WARNING CASES	IF DOOR IS OPEN DURING PROCESS	IF DOOR IS OPEN DURING PROCESS