






- AIR STERILIZATION DEVICE
- H<sub>2</sub>O<sub>2</sub> PLASMA STERILIZATION DEVICE
- HDROZONE PLASMA STERILIZATION DEVICE
- HOSPITAL TYPE EO STERILIZATION DEVICE
- INDUSTRIAL TYPE ETHYLENE OXIDE STERILIZATION DEVICE
- PASSBOX




## PRODUCT ORDERING LIST


**TEKNOMAR**<sup>®</sup>  
Since 1993 **Group**


**We Produce For Your Health...**


PRODUCT	BRAND/ MODEL	DESCRIPTION	SIZE (cm / kg)	PACKAGING SIZE (cm / kg)	ORDER NO
	<b>UVION PATHOGEN IONIZER MEDICAL DEVICE TYPE AIR STERILIZER</b> UVION HOME TYPE	Effective Solution to Eliminate All Airborne Viruses  95% ≤ Bio Decontamination within the First Hour  200 m <sup>3</sup> / h capacity	630 x 150 x 325 cm 10 kg	680 x 200 x 380 cm 11,5 kg	<b>35152</b>
	<b>UVION PATHOGEN IONIZER MEDICAL DEVICE TYPE AIR STERILIZER</b> UVION PROFESSIONAL TYPE	Effective Solution to Eliminate All Airborne Viruses  97% ≤ Bio Decontamination within the First Hour  220 m <sup>3</sup> / h capacity	630 x 150 x 325 cm 10 kg	680 x 200 x 380 cm 11,5 kg	<b>35091</b>
	<b>XRETON-P UVION MOBILE AIR STERILIZER</b>	Atmospheric Plasma Pathogen Ionizer, HEPA Filtered, Mobile Air Sterilizer  96% ≤ BioDecontamination within the First Hour  400 m <sup>3</sup> /h capacity	320 x 320 x 820 cm 28 kg	350 x 350 x 870 cm 29,5 kg	<b>35107</b>




PRODUCT	BRAND/ MODEL	DESCRIPTION	SIZE (cm / kg)	PACKAGING SIZE (cm / kg)	ORDER NO
	<b>UVION AIR ASEPTIZER CW650 IO</b>  <b>ISOLATED ROOM TYPE</b>	Positive-Negative Pressure Air Conditioning Air Purifier Particle Reduction And Decontamination Device BTU: <b>10000</b> Total air change capacity <b>1200 m³/h</b>	79x 58x 229 cm 147 kg	86 x 65 x 237 cm 167 kg Pallet: 92x70 cm	<b>35145</b>
	<b>UVION AIR ASEPTIZER CW650</b>  <b>Standart</b>	Positive-Negative Pressure Air Conditioning Air Purifier Particle Reduction And Decontamination Device BTU: <b>10000</b> Total air change capacity <b>1200 m³/h</b>	79x 58x 229 cm 147 kg	86 x 65 x 237 cm 167 kg Pallet: 92x70 cm	<b>35138</b>
	<b>UVION AIR ASEPTIZER CW1000</b>	Positive-Negative Pressure Air Conditioning Air Purifier Particle Reduction And Decontamination Device BTU: <b>24000</b> Total air change capacity <b>2000 m³/h</b>	109x58x229cm 190 kg	117 x 65 x 237 cm 210 kg Pallet: 122 x 70 cm	<b>35114</b>
	<b>UVION AIR ASEPTIZER CW1350</b>	Positive-Negative Pressure Air Conditioning Air Purifier Particle Reduction And Decontamination Device BTU: <b>24000</b> Total air change capacity <b>2500 m³/h</b>	139 x 58 x 229 cm 217 kg	147 x 65 x 237 cm 237 kg Pallet: 152 x 70 cm	<b>35121</b>
	<b>UVION AIR ASEPTIZER CW1000 CONTAINER</b>	Box Type Positive-Negative Pressure Air Conditioning Air Purifier Particle Reduction And Decontamination Device BTU: <b>24000</b> Total air change capacity <b>2000 m³/h</b>	80 x 139 x 120 cm 200 kg	90 x 150 x 135 cm 230 kg Pallet: 100 x 160 cm	<b>35220</b>

	<b>SAHARAN ASEPTIZER – XRETON-P CW 5000 SAHARAN OUT DOOR</b>	Saharan Type Mobile Air Conditioning and Air Sterilization Device BTU: <b>62000 / 56000</b> Total air change capacity <b>5000 m³/h</b>	334 x 150 x 140 cm 1300 kg		<b>35237</b>
	<b>SAHARAN ASEPTIZER – XRETON-P CW 10000 SAHARAN</b>	Saharan Type Mobile Air Conditioning and Air Sterilization Device BTU: <b>124000 / 112000</b> Total air change capacity <b>10000 m³/h</b>	334 x 150 x 140 cm 1520 kg		<b>35244</b>
	<b>PASSBOX</b>	<b>PASSIVE (STATIC)</b> UV: OPTIONAL INTERLOCK  PassBox in: 60x60x60 cm	64,5 x 84 x 66 cm 48 Kg	75 x 95 x 80 cm 65 Kg Pallet: 120 x 80 cm	<b>35282</b>
	<b>PASSBOX</b>	<b>ACTIVE (DYNAMIC)</b> UV: OPTIONAL INTERLOCK HEPA FILTER  PassBox in: 60x60x60 cm	86 x 87 x 107 cm 98 Kg	96 x 100 x 120 cm 115 Kg Pallet: 120 x 80 cm	<b>35299</b>

PRODUCT	BRAND / MODEL	DESCRIPTION	SIZE (cm / kg)	PACKAGING SIZE (cm / kg)	ORDER NO
	<b>ETO - C 1445 ETHYLENE OXIDE STERILIZATION DEVICES</b>	Hospital Type Ethylene Oxide Gas Sterilizer 134 Lt	98 x 103 x 82 cm 208 kg	110 x 110 x 102 cm 228 kg Pallet: 115 x 115 cm	<b>36050</b>

PRODUCT	BRAND/ MODEL	DESCRIPTION	SIZE (cm / kg)	PACKAGING SIZE (cm / kg)	ORDER NO
	<b>HRF 3000 HYDROGEN PEROXIDE H<sub>2</sub>O<sub>2</sub> PLASMA STERILIZATION DEVICE S75</b>	Single Door Chamber: 113 Lt	73x101x186 cm 355 kg	85 x 110x 195 cm 375 kg Pallet: 90x115 cm	<b>30119</b>
	<b>HRF 3000 HYDROGEN PEROXIDE H<sub>2</sub>O<sub>2</sub> PLASMA STERILIZATION DEVICE S125</b>	Single Door Chamber: 154 Lt	73x101x186 cm 365 kg	85x110x195 cm 385 kg Pallet: 90x115 cm	<b>30133</b>
	<b>HRF 3000 HYDROGEN PEROXIDE H<sub>2</sub>O<sub>2</sub> PLASMA STERILIZATION DEVICE D125</b>	Double Door Chamber: 154 Lt	77x93x185 cm 405 kg	85x100x195 cm 425 kg Pallet: 90x105 cm	<b>30140</b>
	<b>HRF 3000 HYDROGEN PEROXIDE H<sub>2</sub>O<sub>2</sub> PLASMA STERILIZATION DEVICE S175</b>	Single Door Chamber: 192 Lt	82 x 95,5 x 186 cm 450 kg	95 x 101 x 208 cm 520 kg Pallet: 90x115 cm	<b>30157</b>

PRODUCT	BRAND/ MODEL	DESCRIPTION	SIZE (cm / kg)	PACKAGING SIZE (cm / kg)	ORDER NO
	<b>HRF 3000 HDROZONE PLASMA STERILIZATION DEVICE</b> 3S125 LUMEN: 5 M	<b>HDROZONE</b> Hydronium Heavy Molecule & H <sub>2</sub> O <sub>2</sub> Plasma Sterilization Device  <b>Sterilization lumen length capacity: Ø 2 mm x 5 m</b>	Chamber: 154 LT 73 x 95,5 x 198 cm 450 kg	86 x 101 x 208 cm 520 kg Pallet: 90 x 115 cm	<b>30386</b>
	<b>HRF 3000 HDROZONE PLASMA STERILIZATION DEVICE</b> 3S125 LUMEN: 7,5 M	<b>HDROZONE</b> Hydronium Heavy Molecule & H <sub>2</sub> O <sub>2</sub> Plasma Sterilization Device  <b>Sterilization lumen length capacity: Ø 2 mm x 7,5 m</b>	Chamber: 154 LT 73 x 95,5 x 198 cm 450 kg	86 x 101 x 208 cm 520 kg Pallet: 90 x 115 cm	<b>35305</b>
	<b>HRF 3000 HDROZONE PLASMA STERILIZATION DEVICE</b> 3S125 LUMEN: 10 M	<b>HDROZONE</b> Hydronium Heavy Molecule & H <sub>2</sub> O <sub>2</sub> Plasma Sterilization Device  <b>Sterilization lumen length capacity: Ø 2 mm x 10 m</b>	Chamber: 154 LT 73 x 95,5 x 198 cm 450 kg	86 x 101 x 208 cm 520 kg Pallet: 90x115 cm	<b>35312</b>
	<b>HRF 3000 HDROZONE PLASMA STERILIZATION DEVICE</b> 3S125 LUMEN: 15 M	<b>HDROZONE</b> Hydronium Heavy Molecule & H <sub>2</sub> O <sub>2</sub> Plasma Sterilization Device  <b>Sterilization lumen length capacity: Ø 2 mm x 15 m</b>	Chamber: 154 LT 73 x 95,5 x 198 cm 450 kg	86 x 101 x 208 cm 520 kg Pallet: 90x115 cm	<b>35329</b>
	<b>HRF 3000 HDROZONE PLASMA STERILIZATION DEVICE</b> 3S125 LUMEN: 50 M	<b>HDROZONE 50L</b> Hydronium Heavy Molecule & H <sub>2</sub> O <sub>2</sub> Plasma Sterilization Device  <b>Sterilization lumen length capacity : Ø 2 mm x 50 m</b>	Chamber: 154 LT 73 x 95,5 x 198 cm 450 kg	86 x 101 x 208 cm 520 kg Pallet: 90x115 cm	<b>35336</b>

PRODUCT	BRAND/ MODEL	DESCRIPTION	SIZE (cm / kg)	PACKAGING SIZE (cm / kg)	ORDER NO
	<b>HYDROGEN PEROXIDE (H<sub>2</sub>O<sub>2</sub>) SOLUTION CARTRIDGE</b>	HRF 3000 A150 H <sub>2</sub> O <sub>2</sub> Cartridge Sterilizing Agent Hydrogen Peroxide Sterilizing Agent 150 cc	FULL / EMPTY Min. order quantity: 15 pcs in 1 box	30 x 22 x 20,5 cm 4 Kg	<b>30096</b>
	<b>H<sub>2</sub>O<sub>2</sub> C1080 Chemical Indicator For H<sub>2</sub>O<sub>2</sub> &amp; hdrOzone</b>	Hydrogen Peroxide Plasma Sterilizer & HdrOzone Device Chemical Indicator	250 pcs / 1 box	Packing: 17 x 17 cm 250 gr	<b>30000</b>
	<b>EO CARTRIDGE</b>	100% EO GAS	100 gr / 1 cycle		<b>30001</b>



## **TERMS & CONDITIONS**

- I. Prices are given on Euro bss as Ex-works / Ankara deliver term in terms of INCOTERMS 2020.
- II. If the product is available in the stock, delivery duration is 7 business days after order placement; In case of non-availability stock of the product, delivery duration is 4-6 weeks after order placement.
- III. Letter of Credit payment terms are accepted only for order placements above 100.000 € EURO. Below this amount, T/T terms are accepted.
- IV. Devices are under warranty of ONE year against factory faults excluding user faults and consumables.
- V. Spare parts and consumables required for the devices are guaranteed to be supplied for a period of 10 years, subject to a fee
- VI. Teknomar Company reserves rights to make any changes in the field of technological developments.

**TEKNOMAR®**  
Since 1993  
Group



**WE PRODUCE FOR YOUR HEALTH!**

+90 (312) 385 0040

[www.teknomar.com.tr](http://www.teknomar.com.tr)

[info@teknomar.com.tr](mailto:info@teknomar.com.tr)