



Hydrogen Peroxide (H₂O₂) Solution Cartridge

Datasheet

PRODUCT NAME	Hydrogen Peroxide Solution Cartridge
MODEL	HRF 3000 A150
CLASS	IIb
VOLUME	150 ml
SHELF LIFE	6 Months
COMPATIBILITY	HRF 3000 Hydrogen Peroxide Plasma Sterilization Devices HRF 3000 Hydrogen Peroxide Ozone Plasma Sterilization Devices
RFID	Included
CERTIFICATES	ISO 9001.2015 – ISO 13485.2016 – EN ISO 14937 – MDD 93/42/EC-CE Certificate
MATERIAL	PVA/ PVC
INGREDIENT/CHEMICALSUBSTANCE INFORMATION	
FORMULA	Hydrogen Peroxide (H ₂ O ₂)
IDENTITY	(CAS Number) 7722-84-1 (EC Number) 231-765-0 (EC Index Number) 008-003-00 -9 (IUPAC Number) Hydrogen Peroxide
CONCENTRATION	%35-60
DENSITY	1,45 g/cm ³
MOLAR MASS	34,0147 g/mol
MELTING POINT	-0,43 °C
BOILING POINT	150,2 °C